# ODOFIN™

*TasteStrips*


Taste Test 50

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Deutsch .....	3
English .....	10
Symbolerklärung / Definition of symbols .....	17
Auswertung / Evaluation .....	19

PZN: 17975303

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Uetersener Straße 6

D - 25488 Holm

+49 4103 80076-0

+49 4103 80076-29

info@burghart-mt.de

www.burghart-mt.de

**Deutsch**

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**Literatur**

C Mueller, S Kallert, B Renner, K Stiasny, A F P Temmel, T Hummel, G Kobal

**Quantitative assessment of gustatory function in a clinical context using impregnated "taste strips"**

Rhinology. 2003 Mar;41(1):2-6.

PMID: 12677732

Mueller CA, Pintscher K, Renner B.

**Clinical test of gustatory function including umami taste**

Ann Otol Rhinol Laryngol. 2011.

PMID: 21774441

Landis BN, Welge-Luessen A, Brämerson A, Bende M, Mueller CA, Nordin S, Hummel T.

**"Taste Strips" - a rapid, lateralized, gustatory bedside identification test based on impregnated filter papers.**

J Neurol. 2009 Feb;256(2):242-8.

PMID: 19221845.

**English**

**Content**

Contraindication ..... 12

Information for patients ..... 12

Description of the test ..... 12

Labeling of containers and strips ..... 12

Information for examiners ..... 12

Test conditions and preparations ..... 13

Presentation of strips ..... 13

Test sequence ..... 13

Whole-mouth testing ..... 13

Side-separated testing ..... 13

Patient's evaluation ..... 14

Evaluation ..... 14

Shelf life ..... 15

Storage ..... 15

Disposal ..... 15

Accessoires ..... 15

Substances used ..... 15

Literature ..... 16

Symbolerklärung / Definition of symbols ..... 17

Auswertung/Evaluation ..... 19

**IMPORTANT  
 READ CAREFULLY BEFORE USE  
 KEEP FOR FUTURE USE**

### Intended use

The "Taste Strips" test shall be used for psychophysical examination of the sense of taste.

It is meant for detection and quantification of taste disorders.

### User

The test may only be performed by healthcare professionals.

### Age group

There are no age restrictions for the test.

### Warnings

Ask your patients if they are pregnant. The taste quality bitter must NOT be used on pregnant women. Quinine is used as the bitter substance.

Use the tasting strip only ONE time per patient or tongue area application. Afterwards, the taste strip must be disposed of.

If you are allergic to any of the substances used (see Substances used page 15), do not use this test.

Avoid skin contact during administration and use the enclosed tweezers.

Avoid rapid movements with the Taste Strip during application in the oral cavity.

Note that a compromised product can lead to incorrect diagnosis.

The product can be compromised by:

- Change in concentration caused by incorrect storage
- Incomplete test or mix-up of the Taste Strips if the storage containers supplied are not used
- Change in concentration due to use beyond the expiry date
- Failure to apply correctly as described in the instructions for use
- If a Taste Strip is used twice (side-separated testing)
- Incorrect application due to illegible or missing labelling

## Contraindication

Do not use the bitter "Taste Strips" if the patient is pregnant.  
Do not use in case of allergy to any of the substances used.

## Information for patients

Please inform your patient in the following manner: Your sense of taste will be tested. During this test you will be asked to identify the taste of impregnated paper strips. The strips are impregnated with sweet, sour, salty or bitter taste (optionally with umami). Some strips may contain no taste. The substances are not harmful in their used concentrations. Comparable concentrations can be found in food or drinks. If you cannot taste anything, please indicate this ("no taste"). The evaluation of the test takes into account whether you can detect or identify a quality. The test takes about 20 minutes per patient.

## Description of the test

The "Taste Strips" test is a validated examination method to determine the ability to taste. For the assessment of whole-mouth gustatory function the strips are placed on the tongue and the mouth should be closed (the tongue may be moved). In case of regional testing the strips stay in contact with the tongue while the mouth is open until the patient's response.

## Labeling of containers and strips

The container's cap and the strips are labeled with letters A to P and U + V:

highest concentration	A	E	I	M	U + V	Q
	B	F	J	N		R
	C	G	K	O		S
lowest concentration	D	H	L	P		T
	sweet	sour	salty	bitter	blanks	umami

(optional)

## Information for examiners

Some strips of the test are not impregnated with taste solutions. In case of absent sensation, the patient has to answer with "no taste". This procedure is required to document correct identification as well as the ability to detect each quality.

### Test conditions and preparations

Ensure that the patients only drink water, do not smoke and do not eat or use bubble gum etc. one hour before testing. Put the test in front of you on the table. Put a glass of water (low in salt/minerals, not carbonated) and a waste bin in front of the patient. Open the box and take out the container with strips you like to present.

### Presentation of strips

Take the tweezers in one hand and open the cap of the container with the other one. Use the tweezers to pull out one taste strip. In case of touching the strips in the container with bare fingers, the whole container needs to be regarded as contaminated and must be thrown away. The strip is single use only.

### Test sequence

The different taste strips are presented successively in a pseudorandomized manner. After each placement of a strip, the mouth should be rinsed with water. In addition to the impregnated strips, strips with no taste ("Blanks U+V") can be offered. As a recommendation, while the first strip is presented, the patient is told to observe the taste of the paper.

Test sequence I: U-D-P-L-H-U-G-C-O-K-V-F-B-J-N-A-E-I-M

Test sequence II: U-H-L-P-D-U-K-O-C-G-V-N-B-F-J-E-M-I-A

With umami (optional):

Test sequence I: U-D-P-L-T-H-U-G-C-S-O-K-V-F-B-J-N-R-A-E-I-M-Q

Test sequence II: U-H-L-P-T-D-U-S-K-O-C-G-V-N-B-F-R-J-E-Q-M-I-A

### Whole-mouth testing

Place the taste strip in the patient's hand with the following instruction: "Please place the strip in the middle of the tongue, close your mouth and move your tongue and the strip slowly back and forth so that the substances dissolve in the saliva. The strip can also be flipped and placed on the tongue again." After testing a quality, the mouth must be rinsed with water.

### Side-separated testing

In side-separated testing, the right and left sides of the tongue are tested independently of each other. The strips are used like with "whole-mouth testing", but in this case the investigator applies and moves a taste strip on each area to the right or left of the center line of the tongue. The mouth remains open. After testing one side and quality, the mouth must be rinsed with water.

**ATTENTION:** For each side of the tongue a NEW strip (same quality) must be used.

### Patient's evaluation

5 to 10 seconds after each placement of a strip, ask what taste the patient has perceived. The patient may choose between “sweet”, “sour”, “salty”, “bitter” and “no taste”. Mark the answer on the evaluation sheet (see page (19/20)).

### Evaluation

Each correct answer yields one point. If all 16 strips are identified correctly the maximum score is 16 (4 correct for each taste quality). Blanks are not counted during evaluation. Detected strips with incorrect identification may be used to diagnose patients with parageusia. Patients with severe hypogeusia tend to confound strip qualities even at higher concentrations. Therefore, patients with a clear perception of the highest concentrations (A, E, I, M) cannot be diagnosed as ageusic for the quality in question. Patients with hypogeusia (identification <9) may exhibit ageusia for certain taste qualities.

Normative values for each taste quality are shown in the following table:

Number of correctly identified taste qualities				
Sweet	Sour	Salty	Bitter	Total score
≥2	≥2	≥2	≥1	≥9 (Normogeusia)
1 or false identification	1 or false identification	1 or false identification	false identification	<9 (Hypogeusia)
No sensation	No sensation	No sensation	No sensation	No sensation (complete ageusia)

At a score below the threshold value of 9\*, hypogeusia can be assumed regardless of age. A complete ageusia can be assumed if the highest concentrations of all four taste qualities (A, E, I, M) are not detected.

As a test result, the total taste performance and the performance for each taste quality can be evaluated separately. For example: Hypogeusia (total score = 7) with partial ageusia for sour ("no perception of sour").

\*The threshold value of 9 is derived from the 10th percentile from the norm data study.

The evaluation does not take into account the responses for umami, but only sweet, sour, salty and bitter, as described in the instructions.

The data basis is too small to make qualified statements when the responses of umami are added. The additional flavour of umami is therefore only used as a reference so that this area of taste perception is considered but is not included in the evaluation.

### Shelf life

The inside of the cardboard lid, which contains the tasting strips cans, and each tasting strip is marked with a best before date. The format is DD MM YY. Unopened, the tasting strips have a shelf life of 5 years from production. Once opened, the shelf life is limited to 12 months. Enter the period of the first opening, on the symbol located outside the box lid.

### Storage

Keep dry, at room temperature (10 – 30 °C), do not put into fridge.

### Disposal

Used strips must be disposed of immediately in the household waste.

### Accessoires

Component	Item number
Tweezers	LA-13-00310
Storage container	LA-13-00445
Instructions for use and evaluation material	LA-18-00051

### Substances used

Taste quality	Substance
sour	Citric acid
salty	Sodium chloride
sweet	Sucrose
bitter	Quinine hydrochloride dihydrate
Umami (optional)	Sodium L-glutamate monohydrate
Solvent	Aqua Bidest

**Literature**

C Mueller, S Kallert, B Renner, K Stiassny, A F P Temmel, T Hummel, G Kobal

**Quantitative assessment of gustatory function in a clinical context using impregnated "taste strips"**

Rhinology. 2003 Mar;41(1):2-6.

PMID: 12677732

Mueller CA, Pintscher K, Renner B.

**Clinical test of gustatory function including umami taste**

Ann Otol Rhinol Laryngol. 2011.

PMID: 21774441

Landis BN, Welge-Luessen A, Brämerson A, Bende M, Mueller CA, Nordin S, Hummel T.



**"Taste Strips" - a rapid, lateralized, gustatory bedside identification test based on impregnated filter papers.**

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