ODOFIN™

Sniffin' Sticks

Discrimination Test

burgardt
MESSTECHNIK

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ENGLISH

Content

Description smell test with odour sticks 12
Intended use 12
User 12
Age group 12
Warnings 12
Contraindication 13
Patient's information for the smell test 13
Usage of the Sniffin' Stick 14
Execution Discrimination Test 14
Evaluation of the results 15
Shelf life 15
Storage 15
Disposal 15
Ingredients 16
Example 16
Literature 17
Evaluation discrimination test (copy template) 29
Normative data 30

**IMPORTANT
READ CAREFULLY BEFORE USE
KEEP FOR FUTURE USE**

Description smell test with odour sticks

ODOFIN™ „Sniffin` Sticks“ are used to investigate human olfactory performance. They have been developed in close cooperation with the association „Olfaktologie und Gustologie“ of the German Society for Ear, Nose and Throat Medicine, Head and Neck Surgery.

Intended use

Discrimination of odours is based on comparison of 3 smell presentations (triplet). To this end, two times the same odour is presented (non-target) and one time a different smelling one (target). The subject's task is to indicate which one smells different. This comparison is performed for 16 triplets.

Name	Font colour	Colour of stick's end	Number of sticks
Discrimination Test	Green	Green (target) Red, blue (non-target)	16 Triplets = 48 sticks

User

The test is to be performed by healthcare professionals only.

Age group

The test can be carried out starting at the age of 5 years.

Warnings

The olfactory sticks must not touch the patient's skin during the olfactory presentation. However, if the stick is inadvertently touched, it must be disposed of for hygienic reasons and because of the reduced shelf life of the odours and replaced by a new one.

All substances used in this test are not health hazardous in their used concentrations. Rarely allergic reactions show up after inhaling and skin contact. In case of allergic reaction contact the attending doctor and the manufacturer.

The user shall wear odour-free gloves during the procedure.

Use of the sticks only in an odourless or well-ventilated environment.

Only use the stick for recommended procedures and not for writing purposes.

Caution: Ask for intolerances before each test.

Store the product safely and out of the reach of patients and children.

If the tip of the Sniffin' Sticks accidentally comes in contact with the gloves, the gloves need to be replaced to avoid contamination of other odours.

Contraindication

An allergy against one of the odorants used in the tests.

Patient's information for the smell test

The patient shall be informed about the test in the following way:

Your sense of smell is to be tested. This will be done by means of odour sticks "Sniffin' Sticks". You will be blindfolded and smell the sticks which contain different odorants. The test will take approximately 10-15 minutes. All substances contained in the test are non-toxic and not harmful in the used concentrations. If you should not be able to smell and thus cannot decide on an answer, please just guess. The evaluation of the test will take into account that you have been asked to make a decision in any case.

Question: How well can you distinguish smells from one another? Repeatedly, three sticks will be presented to you one after another. Two of these contain the same odorant; a third one contains a different smell. You have to indicate the stick that smells differently.

Information for the examiner:

Patients often do not understand why they have to decide for a particular Sniffin' Stick, even though they did probably not smell anything. The patient may be given the additional explanation that this is a basic principle of modern test procedures and that accidentally correct answers are taken into account during evaluation of the test. During the evaluation it will be taken into account that they usually reach a score greater than zero, even though they did not actually smell anything, but just hit the correct answer by chance. Thus, reaching this score indicates that the patient cannot smell or cannot properly smell. The Discrimination test may only be used in combination with instructions for use, holders and sticks. Only the refill set supplied by the manufacturer may be used as a refill set.

Usage of the Sniffin' Stick

When testing, the cap of only one stick should be removed. The cap has to be put back on the odour stick right after presenting the odour. Afterwards, the odour stick should immediately be placed and turned into the holder, with the tip pointing downwards.

For testing, the opened odour stick is held 2 cm centred in front of both nostrils and the patient is asked to smell e.g. using the sentence "Please sniff now". Each stick should not be presented for more than 3-4 seconds. The interval between the different presentations should be around 30 seconds.

If the test is to be carried out left or right sided, the opened stick is simply held in front of one nostril. The untested side is in this case closed by the patient by placing the pad of the right or left thumb on the nostril from below. The closure must not deform the nose.

To gain results as significant as possible the patients should not have anything but water for at least 30 minutes before the test. They should also refrain from consuming chewing gum, sweets or cigarettes. A short anamnesis can be recorded in the evaluation form. This includes: Name, age, occupation and smoking habits of the patient, as well as the information whether the test was made bilateral (L+R), left sided (L) or right sided (R). The patients should not receive any hint of how correct their answers are during the test.

Execution Discrimination Test

The patients are blindfolded by putting on an "eye cover". 16 triplets (with three sticks each) are presented to the patients. Two sticks of the triplet smell the same (non-target), only one stick of this triplet contains a different odorant (target). It is the one with the green cap. It is the patient's task to identify this different smelling stick (target). **The patient always has to make a choice, even if not sure.** The interval between the first stick of one triplet and the first stick of the next triplet should be 30 seconds. Each stick is presented only once. Repeated presentation of a triplet on demand of the patient is not admissible. The examiner always holds three sticks while testing. In the process only the three sticks should be taken from the holder, which are currently needed. The presentation order of the stick

Presentation	sequence of the sticks
1	red – green – blue
2	blue – red – green
3	green – blue – red
4	red – green – blue
Repeat the pattern until the test is over!	

containing the odorant has to be alternated by the examiner. This cycle is repeated throughout the whole test.

1. The colour of the stick named by the patient as the one smelling different (target) is marked at the corresponding spot in the enclosed evaluation sheet (see page 29).
2. At the end of the test, the number of correctly detected triplets is added to obtain the score.

The procedure is explained again under example (see page 16).

Evaluation of the results

Since the Discrimination Test by itself is not intended for a detailed assessment, the only possible statement that can be made is based on the norm values on page 30. If the result is within the range of the norm values, there is no hint to a smell disorder of the patient.

The result of the Discrimination Test is usually evaluated within the scope of the Extended Test together with the Threshold and the Identification Test.

Shelf life

The sticks are labelled with a lot-number and a best before date. Use beyond the expiry date may lead to incorrect results.

Storage

Refer to the symbols for how to store the test (dry, protect from sunlight). The temperature should be between 10°C and 30°C.

The sticks must be stored in an upright "overhead" position (white cap downwards).

Storage in the fridge can lead to a change in the intensity of the odour.

Disposal

Unusable sticks must be disposed of in their sealed state in the household waste.

Ingredients

Stick	Chemical Target	Chemical non Target (2x)
No. 1	n-Octyl acetate	Cinnamaldehyde
No. 2	n-Butanol	2-Phenylethanol
No. 3	Isoamyl acetate	Anethole
No. 4	Anethole	Eugenol
No. 5	Geraniol	n-Octyl acetate
No. 6	2-Phenylethanol	Isoamyl acetate
No. 7	(+)-Limonene	(+)-Fenchone
No. 8	(-)-Carvone	(+)-Carvone
No. 9	(-)-Limonene	Citronellal
No. 10	2-Phenylethanol	L-Menthol
No. 11	(+)-Carvone	Geraniol
No. 12	n-Butanol	(-)-Fenchone
No. 13	Citronellal	Linalool
No. 14	Pyridine	(-)-Limonene
No. 15	Eugenol	Cinnamaldehyde
No. 16	Eucalyptol	α-Ionon
Solvent	Propylenglycol	

Example

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Red	x									x						x
Green		x	x		x	x			x		x	x	x	x	x	
Blue				x			x	x								

The green stick is the target, i.e. this stick must be detected because it is different from the other two sticks.

In this example, the patient detected the target 10 times and reached a score of 10.

Literature

A Oleszkiewicz , V A Schriever , I Croy , A Hähner , Thomas Hummel

Updated Sniffin' Sticks normative data based on an extended sample of 9139 subjects

Eur Arch Otorhinolaryngol. 2019 Mar;276(3):719-728. doi: 10.1007/s00405-018-5248-1. Epub 2018 Dec 15.

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IMPORTANT
READ CAREFULLY BEFORE USE
KEEP FOR FUTURE USE

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
Evaluation discrimination test (copy template)

Date, Time _____, _____

Name _____ Age _____

Occupation _____

Smoker Yes No

Discrimination test

left-sided testing

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Red																
Green																
Blue																

right-sided testing

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Red																
Green																
Blue																

bilateral testing

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Red																
Green																
Blue																

Result (Sum of the correct Discriminations*)

left _____ right _____ bilateral _____

* Green is correct

Normative data


The following graph can be used for the evaluation. The graph was created with the data from the norm data study mentioned under "Literature".

The classification of the ranges normosmia and hyposmia, results from the data of the most powerful age groups (21-30) with regard to the ability to smell.

Since the boundaries for anosmia and hyposmia coincide according to the results of the validation studies, this test alone cannot be used to differentiate between anosmia and hyposmia.

The graph is not separated by sex, as the difference of the olfactory performance between the sexes is not of a relevant magnitude for the evaluation.

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