ODOFIN™

Sniffin' Sticks

Extended Test

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English

Content

Intended use	15
User	15
Age group	15
Warnings	15
Contraindication	15
Patient information on the smell test	16
Overview and labelling of the subtests	17
Usage	17
Application of the Threshold Test.....	18
Application of the Discrimination Test	19
Identification Test.....	19
Evaluation of the Results	19
Shelf Life.....	20
Storage	20
Disposal	20
Ingredients	21
Example.....	22
Literature	24
Evaluation Extended Test	42
TDI – Score.....	43
Evaluation Threshold Test	44
Evaluation Discrimination Test.....	45
Evaluation Identification Test	46

IMPORTANT
READ CAREFULLY BEFORE USE
KEEP FOR FUTURE USE

Intended use

The Extended Test (also called "TDI Test" - Threshold, Discrimination, Identification) consists of 3 individual tests. Each of these tests generates a score, all of are added up to the TDI score at the end. Based on this score, the user can determine a diagnostic classification of the olfactory ability.

User

The test is to be performed by healthcare professionals only.

Age group

The test can be carried out starting at the age of 5 years.

Warnings

The olfactory sticks must not touch the patient's skin during the olfactory presentation. However, if the stick is inadvertently touched, it must be disposed of for hygienic reasons and because of the reduced shelf life of the odours and replaced by a new one.

All substances used in this test are not health hazardous in their used concentrations. Rarely allergic reactions show up after inhaling and skin contact. In case of allergic reaction contact the attending doctor and the manufacturer.

The user shall wear odour-free gloves during the procedure.

Use of the sticks only in an odourless or well-ventilated environment.

Only use the stick for recommended procedures and not for writing purposes.

Caution: Ask for intolerances before each test.

Store the product safely and out of the reach of children.

If the tip of the Sniffin' Sticks accidentally comes in contact with the gloves, the gloves need to be replaced to avoid contamination of other odors.

Contraindication

An allergy against one of the odorants used in the tests.

Patient information on the smell test

The patient should be informed about the test in the following way:

Your sense of smell should be tested. This is done with the help of smelling sticks "ODOFIN™ Sniffin' Sticks". You smell sticks containing different odour components sometimes with open eyes and sometimes blindfolded. The test is divided into three parts, it takes about 30-40 minutes. All substances used in the test are non-toxic and non-hazardous in the used concentrations. If you can't smell anything and can't decide which answer to give, just guess. When evaluating the test, it is taken into account that you were asked to make a decision in each case.

• Threshold Test:

Question: At which concentration can you perceive an odorant? Repeatedly, three sticks will be presented to you one after another. Only one of these sticks contains the odorant, the other two contain a solvent. You have to find out the stick in which you think you have perceived the odorant. In any case, you have to pick one of the three sticks, even if you are not sure of your decision.

• Discrimination Test:

Question: How well can you distinguish smells from one another? Repeatedly, three sticks will be presented to you one after another. Two of these contain the same odorant, a third one contains a different smell. You have to indicate the stick that smells differently. In any case, you have to pick one of the three sticks, even if you are not sure of your decision.

• Identification Test:

Question: How well can you identify smells? 16 smells known from everyday life will be presented to you. You have to assign each smell to one of four terms on a card that describes the smell best.

Information for the Examiner:

Patients often do not understand why they have to decide for a particular smell on the answer card, even though they did probably not smell anything. The patient may be given the additional explanation that this is a basic principle of modern test procedures and that accidentally correct answers are taken into account during evaluation of the test. During the evaluation it will be taken into account that they usually reach a score greater than zero, even though they did not actually smell anything, but just hit the correct answer by chance. Thus, reaching this score indicates that the patient cannot smell or cannot properly smell.

Requirement

- The total duration of the test is approx. 30-40 minutes.
- To obtain the most reliable results possible, the patient should consume nothing other than water at least 30 minutes before the start of the test.
- Chewing gum, sweets or cigarettes must also be avoided.
- During the test, the patient is not to be given any indication of the correctness of the decision.
- The Extended Test may only be used in combination with instructions for use, holders and sticks. Only the refill set supplied by the manufacturer may be used as a refill set.

Overview and labelling of the subtests

Threshold Test

The threshold test is intended to determine the olfactory threshold of a subject with the help of a graduated concentration of n-butanol. A forced-choice staircase procedure is used, whereby the patient must make a statement which determines a score with the help of turning points.

Discrimination Test

The differentiation of odours is based on a comparison between 3 odours (triplets). Therefore the same odour is offered twice (non-target) and another odour (target). The patient must indicate the odour that smells different in each case. These comparisons are performed for 16 triplets.

Identification Test

Here the ability to identify everyday smells by means of a card with 4 choices is determined. This is a "multiple-forced-choice" procedure in which the patient must select one of the 4 terms. During the identification test a total of 16 odours are presented.

The different subtests of the Extended Test are marked as follows:

Name	Font Colour	Colour of Stick's End	Number of Sticks
Threshold Test	Red resp. Purple	Red (target) Green, blue (blanks)	16 triplet = 48 sticks
Discrimination Test	Green	Green (target) Red, blue (non-target)	16 triplet = 48 sticks
Identification Test	Black	Blue	16 sticks

Usage of the Sniffin' Stick

The 3 parts of the Extended Test should be applied in the following order:

1. Threshold Test
2. Discrimination Test and
3. Identification Test.

A break of about three to five minutes should be taken between the three tests.

During the examination of the olfactory threshold and the discrimination ability, the patient must be blindfolded. A surgical mask or similar can be used for this purpose.

For odour testing, only the cap of one stick has to be removed. The cap must be put back on the odour stick immediately after the odour presentation. The odour stick must then be immediately placed back into the holder with the cap facing down. The noises during the test should be independent of the patient's answer (i.e., even if the answer is correct during the threshold test, the sticks must be inserted back into the holder to give no acoustic indication).

For testing, the opened odour stick is held 2 cm centred in front of both nostrils and the patient is asked to smell, e.g. using the sentence "Please sniff now". Each stick should be presented for 3-4 seconds. The interval between the different presentations should be around 30 seconds.

If the test is to be carried out left or right sided, the opened stick is simply held in front of one nostril. The untested side is in this case closed by the patient by placing the pad of the right or left thumb on the nostril from below. The closure must not deform the nose.

A short anamnesis can be recorded in the evaluation form. This includes: Name, age, illnesses, medication, occupation and smoking habits of the patient, as well as the information whether the test was made bilateral (L+R), left sided (L) or right sided (R).

Make sure that the examination is not disturbed.

Application of the Threshold Test

The patient is first familiarised with the smell of the target stick. Therefore the stick with the highest concentration (stick no. 1 with the red stick's end) is presented. For the actual test the patient is blindfolded by putting on an eye cover (e.g., a sleeping mask or surgical mask).

Three sticks (triplet) are each presented to the patient at approximately 5-second intervals. Only one stick of this triplet contains the odorant (the stick with the red cap), the other two contain just the diluent (sticks with blue and green cap = blanks). It is the patient's task to find out the stick that smells different from the two other sticks. If the target stick is recognised correctly, the triplet is presented again (see note below). Each stick is presented only once, repeated presentation of a triplet on demand of the patient is not admissible. The interval between the first stick of one triplet and the first stick of the next triplet should be 30 seconds. The highest dilution level is stick 16 (high dilution = low fragrance concentration, low dilution = high fragrance concentration).

The examiner always holds three sticks while testing. In the process only the three sticks should be taken from the holder, which are currently needed. The presentation order of the stick containing the odorant has to be alternated by the examiner after each triplet. To accomplish this, the three sticks of a triplet are first presented in the order "red", "green", "blue", according to the colour code of the caps. For the next triplet in the order "blue", "red", "green" and then in the order "green", "blue", "red". This cycle is repeated throughout the whole test.

Presentation	sequence of the sticks
1	red – green – blue
2	blue – red – green
3	green – blue – red
4	red – green – blue
The sample repeats itself until the test is over!	

1. When starting the test, triplets are presented in descending order 16, 14, 12, until the patient has recognised a triplet correctly. A certain odorant concentration **is only identified correctly if the stick containing the odorant is recognised twice in a row**, that is, if the stick containing the odorant is identified again in a second presentation of the same triplet (although a triplet is only presented a second time if the patient has identified the odorant correctly during the first presentation). **The patient always has to make a choice, even if not sure - "forced choice"**.
2. This initially correctly identified dilution step is the starting point for the following test procedure. It is marked with two crosses in the first column of the enclosed protocol sheet. After that, the next higher dilution step is presented. This is now documented in the second column.
3. If this is also identified correctly twice, again the next higher dilution step is presented and so on, until the patient makes a wrong decision. Again, the corresponding dilution step, at which the error occurred, is marked with a minus in the protocol sheet, this time in the second column.
4. After that, the next lower dilution step is presented.
5. If this is not identified by the patient, the next lower dilution step is presented and so on, until the patient identifies a dilution step correctly (i.e., twice). Again, this dilution step is marked with two crosses in the protocol sheet. Then a higher dilution step is presented and so on.
6. The test is finished when 7 turning points are passed through. The smell threshold is defined as the average value (result) of the dilution steps which have been marked at the last 4 turning points which were passed through.

The procedure is explained again under Examples (see page 22).

Note: Since the blanks are not empty, but instead contain a diluent, they have a faint characteristic smell. This is intentional, because the dilutions of the target sticks are made with this substance. This way a common background smell exists for target and blanks.

Application of the Discrimination Test

The patient is blindfolded by putting on an eye cover. 16 triplets (with three sticks each) are presented to the patient. Two sticks of the triplet smell the same (non-target), only one stick of this triplet contains a different odorant (target). It is the one with the green cap. It is the patient's task to identify this differently smelling stick (target). **The patient always has to make a choice, even if not sure.** The interval between the first stick of one triplet and the first stick of the next triplet should be 30 seconds. Each stick is presented only once, repeated presentation of a triplet on demand of the patient is not admissible.

The examiner always holds three sticks while testing. In the process only the three sticks should be taken from the holder, which are currently needed. The presentation order of the stick containing the odorant has to be alternated by the examiner within each triplet as in the threshold test (see page 18). The three sticks are each presented to the patient at approximately 5-second intervals.

1. The colour of the stick named by the patient as the one smelling different (target) is marked at the corresponding spot in the enclosed protocol sheet.
2. At the end of the test, the number of correctly detected triplets is added to obtain the score.

The procedure is explained again under Examples (see page 22).

Application of the Identification Test

16 sticks are presented to the patient one after another at an interval of approximately 30 seconds. In order to be able to name the odorants, the patient gets a card with four terms each (multiple-choice card) before the pen is presented. It is the patient's task to pick the term which describes the presented odorant best. **The patient always has to make a choice, even if not sure.** The letter of the term given by the patient is marked in the enclosed protocol sheet (A-D).

A template foil is enclosed with the test for evaluation. Place it on the identification test table and add the correct answers (white field of the template). The sum is the score.

If the test is to be carried out left or right sided, consider that in this case the order of the presented sticks should be altered to reduce the predictability of the answer. It is e.g., possible to start with stick no. 1 at the left nostril and with stick no. 8 at the right nostril. Another option is the alternating application, on the left and right side, until all 16 smells have been presented on both sides. It is advisable to stay with the initially selected procedure and describe it on the enclosed protocol sheet. Whichever method you choose, the important thing is to start identifying on the worse side.

The procedure is explained again under Examples (see page 22).

Evaluation of the Results

The result of this test is expressed as the sum of the results of the 3 subtests, the so called TDI score (threshold, discrimination, identification). Here a score of ≥ 30.75 rates as normal, a score of 30.5 or less indicates hyposmia and a score of 16 and below points to functional anosmia in form of a complete loss of the sense of smell or an extremely weakened smell ability according to the association Olfactology and Gustology of the German ENT-Society.

If no threshold can be determined, i.e., if the patient cannot distinguish the odour in stick 1 from the blanks, the score of the threshold test is set as 1. If the TDI value is ≤ 16 , the patient is functionally anosmic.

In addition, there is a category of supersmellers who achieve an TDI value ≥ 41.5 . These are individuals with an exceptional sense of smell.

The result allows a detailed evaluation of the olfactory function. Since this is a subjective test procedure, objective olfactometry is recommended if fraud or deceit is suspected.

Shelf Life

The sticks are labelled with a lot-number and an expiry date. Use beyond the expiry date may lead to incorrect results.

Printout of the lot-no. and the best before date (example picture)

Storage

Refer to the symbols for how to store the test (dry, protect from sunlight). The temperature should be between 10°C and 30°C.

The sticks must be stored in an upright overhead position (white cap downwards).

Caution:

Storage in the fridge can lead to a change in the intensity of the odour.

Disposal

Unusable sticks must be disposed of in their sealed state in the household waste.

Ingredients

Identification	
Stick	Odour
No. 1	Orange oil
No. 2	Leather fragrance
No. 3	Cinnamaldehyde
No. 4	Peppermint oil
No. 5	Banana fragrance
No. 6	Citrus oil
No. 7	Anethole
No. 8	Turpentine oil
No. 9	Garlic oil
No. 10	Coffee oil
No. 11	Apple flavour
No. 12	Clove oil
No. 13	Pineapple aroma odor
No. 14	Rose perfum oil
No. 15	Anethole
No. 16	Fish flavour
Solvent	Diethyl phthalate
Solvent	Propylenglycolum

Discrimination		
Stick	Chemical Target	Chemical non-Target (2x)
No. 1	n-Octyl acetate	Cinnamaldehyde
No. 2	n-Butanol	2-Phenylethanol
No. 3	Isoamyl acetate	Anethole
No. 4	Anethole	Eugenol
No. 5	Geraniol	n-Octyl acetate
No. 6	2-Phenylethanol	Isoamyl acetate
No. 7	(+)-Limonene	(+)-Fenchone
No. 8	(-)-Carvone	(+)-Carvone
No. 9	(-)-Limonene	Citronellal
No. 10	2-Phenylethanol	L-Menthol
No. 11	(+)-Carvone	Geraniol
No. 12	n-Butanol	(-)-Fenchone
No. 13	Citronellal	Linalool
No. 14	Pyridine	(-)-Limonene
No. 15	Eugenol	Cinnamaldehyde
No. 16	Eucalyptol	α -Ionon
Solvent	Propylenglycolum	

Threshold	
Stick	Chemical
No. 1-16	2-Phenylethanol or n-Butanol
No. 1-16	Propylenglycolum or Aqua conservans

In case of contact with skin, wash thoroughly with water, consult safety data sheets if necessary. The specific safety data sheets can be obtained on request from the manufacturer.

Examples

Threshold test:

X means target detected

- means target not detected

Dilution	Turning points						
	1	2	3	4	5	6	7
1							
2							
3							
4							
5							
6	x x		x x		x x		x x
7	x -	x x	-	-		x x	-
8	x -	x x	x -			x -	
9		x -					
10	-						
11							
12	-						
13							
14	-						
15							
16	-						

1: In this example, the starting point is at the 6th dilution step, since the target stick was detected twice in series. Previously, steps 16, 14, 12, 10 and 8 were not detected or were not detected twice. Step 6 is also the first turning point. (Start with stick 16)

2: The next higher dilution step has now been selected (7). Two presentations were detected, continuing to the next higher dilution step (8). This was repeated until no two target sticks were detected in the 9th dilution step. This means that the second turning point has been found.

3: The next lower dilution step was then selected (8). No two target sticks were detected here, continuing to the next lower dilution level (7). No target stick was detected during the first presentation. In step 6, two target sticks have been detected. This means that the third turning point has been found.

This procedure is repeated until seven turning points have been found.

4: The score is determined from the last four turning points (i.e., turning points 1-3 are not taken into account). The average of the last four turning points is calculated for this purpose. In this example, the patient would have achieved a score of 6.75 $((7+6+8+6)/4=6.75)$.

Note: A turning point always means a change of direction.

Note: If stick 1 is also not recognised during the search for the first turning point, a smelling disorder is assumed and the score is set to 1.

Note: If stick 16 is correctly recognised twice in the first performance, this is the first turning point. As there are no sticks with a higher dilution level, the test is continued with stick 16.

Discrimination test:

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Red	x									x						x
Green		x	x		x	x			x		x	x	x	x	x	
Blue				x			x	x								

The green stick is the target, i.e., this stick must be detected because it is different from the other two sticks.

In this example, the patient detected the target 10 times and reached a score of 10.

Identification test 16:

The thickly-bordered boxes in the table for the left-sided test are intended to represent the patient's answers

Now place the supplied answer template on the identification sheet and add the correct answers. In this example, the patient has reached a score of 10.

ODOFIN™ Sniffin' Sticks

←				
1	A	B	C	D
2	A	B	C	D
3	A	B	C	D
4	A	B	C	D
5	A	B	C	D
6	A	B	C	D
7	A	B	C	D
8	A	B	C	D
9	A	B	C	D
10	A	B	C	D
11	A	B	C	D
12	A	B	C	D
13	A	B	C	D
14	A	B	C	D
15	A	B	C	D
16	A	B	C	D

→				
1	A	B	C	D
2	A	B	C	D
3	A	B	C	D
4	A	B	C	D
5	A	B	C	D
6	A	B	C	D
7	A	B	C	D
8	A	B	C	D
9	A	B	C	D
10	A	B	C	D
11	A	B	C	D
12	A	B	C	D
13	A	B	C	D
14	A	B	C	D
15	A	B	C	D
16	A	B	C	D

↔				
1	A	B	C	D
2	A	B	C	D
3	A	B	C	D
4	A	B	C	D
5	A	B	C	D
6	A	B	C	D
7	A	B	C	D
8	A	B	C	D
9	A	B	C	D
10	A	B	C	D
11	A	B	C	D
12	A	B	C	D
13	A	B	C	D
14	A	B	C	D
15	A	B	C	D
16	A	B	C	D

Σ

10

Σ

Σ

ODOFIN™ Sniffin' Sticks

←				
1	A	B	C	D
2	A	B	C	D
3	A	B	C	D
4	A	B	C	D
5	A	B	C	D
6	A	B	C	D
7	A	B	C	D
8	A	B	C	D
9	A	B	C	D
10	A	B	C	D
11	A	B	C	D
12	A	B	C	D
13	A	B	C	D
14	A	B	C	D
15	A	B	C	D
16	A	B	C	D

→				
1	A	B	C	D
2	A	B	C	D
3	A	B	C	D
4	A	B	C	D
5	A	B	C	D
6	A	B	C	D
7	A	B	C	D
8	A	B	C	D
9	A	B	C	D
10	A	B	C	D
11	A	B	C	D
12	A	B	C	D
13	A	B	C	D
14	A	B	C	D
15	A	B	C	D
16	A	B	C	D

↔				
1	A	B	C	D
2	A	B	C	D
3	A	B	C	D
4	A	B	C	D
5	A	B	C	D
6	A	B	C	D
7	A	B	C	D
8	A	B	C	D
9	A	B	C	D
10	A	B	C	D
11	A	B	C	D
12	A	B	C	D
13	A	B	C	D
14	A	B	C	D
15	A	B	C	D
16	A	B	C	D

Σ

10

Σ

Σ

Literature

A Oleszkiewicz , V A Schriever , I Croy , A Hähner , Thomas Hummel

Updated Sniffin' Sticks normative data based on an extended sample of 9139 subjects

Eur Arch Otorhinolaryngol. 2019 Mar;276(3):719-728. doi: 10.1007/s00405-018-5248-1. Epub 2018 Dec 15.

PMID: 3055435

"ODOFIN™ Sniffin' Sticks" were developed in close cooperation with the working group "Olfactology and Gustology" of the German Society for Otorhinolaryngology, Head and Jaw Surgery.

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Diagnosis, History

Drugs

The patient judges his or her own olfactory sensitivity as

normal

☐

decreased

☐

increased

☐

Results – Extended Test

Threshold

	left
--	------

	right
--	-------

	bilateral
--	-----------

Highest value

--

+

Discrimination

	left
--	------

	right
--	-------

	bilateral
--	-----------

Highest value

--

+

Identification

	left
--	------

	right
--	-------

	bilateral
--	-----------

Highest value

--

=


--

TDI – Score

≥ 41.5	Supersmeller
$41.25 - 30.75$	Normosmia
$30.5 - 16.25$	Hyposmia
≤ 16	Functional anosmia

The following graph can be used for the evaluation. The graph was created with the data from the norm data study mentioned under "Literature". The graph is not separated according to sex, since the difference of the olfactory performance between the sexes is not of a relevant magnitude for the evaluation.

The classification of supersmeller, normosmia and hyposmia is derived from the data of the best-performing age group (21-30) with regard to the ability to smell. The differentiation from anosmia is described in the study.

Example:

If a patient's score is on the 10th percentile, this means that 90% of patients in the same age group have a higher score and 10% a lower one.

Evaluation Threshold Test (copy template)

Date, time _____, _____

Age _____

Name _____

Threshold Test

left-sided testing

Dil.							
1							
2							
3							
4							
5							
6							
7							
8							
9							
10							
11							
12							
13							
14							
15							
16							

right-sided testing

Dil.							
1							
2							
3							
4							
5							
6							
7							
8							
9							
10							
11							
12							
13							
14							
15							
16							

Results (Average of the last four turning points)

left _____ right _____

bilateral testing

Dil.							
1							
2							
3							
4							
5							
6							
7							
8							
9							
10							
11							
12							
13							
14							
15							
16							

Result (Average of the last four turning points)

bilateral _____

Evaluation Discrimination Test (copy template)

Date, time _____, _____

Age _____

Name _____

Discrimination Test

left-sided testing

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Red																
Green																
Blue																

* Green is correct

right-sided testing

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Red																
Green																
Blue																

* Green is correct

Results (Sum of the correct Discriminations*)

left _____ right _____

bilateral testing

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Red																
Green																
Blue																

* Green is correct

Result (Sum of the correct Discriminations*)

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