# ODOFIN™

*TasteStrips*


## Taste Test 25

**burghart**  
MESSTECHNIK

DEUTSCH .....	3
ENGLISH.....	10
DANSK.....	17
LEGENDE / LEGEND / LEGENDE.....	23

PZN: 17975326

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Uetersener Straße 6  
D - 25488 Holm

+49 4103 80076-0

+ 49 4103 80076-29

info@burghart-mt.de

www.burghart-mt.de

CE

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IMPORTANT  
READ CAREFULLY BEFORE USE  
KEEP FOR FUTURE USE

**Intended use**

Orienting examination to see if the ability to taste is present.

**User**

The test may only be performed by healthcare professionals.

**Scope of application**

Approximately 10-15% of patients with olfactory disorders also have taste disorders. It is therefore advisable to check not only the olfactory ability but also the tasting ability. The stripes in the 4 basic taste qualities "sweet", "sour", "salty" and "bitter" allow an examination of the tasting ability.

**The Taste Strips 25 are an accompanying test to the Screening 12 test and can be used as an orientation for further investigations.**

A more detailed examination is possible with Taste Strips 50. This contains the taste strips with basic taste qualities "sweet", "sour", "salty" and "bitter" in 4 different concentrations each, in order to determine the severity of a taste disorder.

**Age group**

There are no age restrictions for the test.

**Warnings**

Ask your patients if they are pregnant. The taste quality bitter must NOT be used on pregnant women. Quinine is used as the bitter substance.

Use the Taste Strip only ONE time per patient or tongue area application. Afterwards, the Taste Strip must be disposed of.

If you are allergic to any of the substances used (see table on page 15), do not use this test.

Avoid skin contact during administration and use the enclosed tweezers.

Avoid rapid movements with the Taste Strip during application in the oral cavity.

Note that a compromised product can lead to an incorrect assessment.

The product can be compromised by:

- Change in concentration caused by incorrect storage.
- Incomplete test or mix-up of the Taste Strips if the storage containers supplied are not used.
- Change in concentration due to use beyond the expiry date.
- Failure to apply correctly as described in the instructions for use.
- Incorrect application due to illegible or missing labelling.

### **Contraindication**

Do not use the bitter "Taste Strip" if the patient is pregnant.  
Do not use in case of allergy to any of the substances used.

### **Patient information**

The patient should be informed about the test in the following way: Your sense of taste will be tested. This is done with the help of taste strips "Taste Strips 25". During this test you will be asked to identify the taste of impregnated paper strips. The strips are impregnated with sweet, sour, salty or bitter taste. The used substances are not harmful in their used concentrations. Comparable concentrations can be found in food or drinks. If you cannot taste anything, please indicate this ("no taste"). The evaluation of the test takes into account whether you can detect or identify a quality.

The test takes about 5 minutes per patient.

### **User information**

The application is done by placing the taste strips on the tongue and closing the mouth (the tongue can be moved by the patient).

If the examiner is interested in the gustatory sensitivity of certain areas of the tongue (e.g. side separated testing), the mouth remains open and the strip is only in contact with this area until the patient can give an answer. The substances used in the test are harmless to health in the presented form and are used in food (sucrose, table salt, citric acid, quinine as e.g. in bitter lemon).

### **Presentation**

Take the tweezers in one hand and open the can with the other. Take out the taste strips **ONLY** with the supplied tweezers.

The patient's answer can be noted on the evaluation sheet (  ) but will NOT affect the score of the Screening 12 test.

## Implementation

Take out the taste strip with the enclosed tweezers and place it in the middle of the tongue so that the taste strip still sticks out a little from the mouth. The patient should now close his mouth and move his tongue slightly. After 5 to 10 seconds the patient can be asked what taste the patient has perceived. Please remove the taste strip again.

It is recommended to wait at least 30 seconds between the presentation of the different strips and to let the patient rinse his mouth with water. Furthermore, "bitter" (can M) should be offered last, as this taste is present for the longest time. As with the Screening 12 test, a selection card with terms is available.

## Shelf life

Each taste strip is provided with a best-before date. The format is DD MM YY. Unopened, the taste strips have a shelf life of 5 years from Production.

After opening, the shelf life is limited to 12 months. Note the date of first opening on the symbol (12M) stuck on the container.

## Storage

Keep dry, at room temperature (10 – 30 °C), do not put into fridge.

## Disposal

Used taste strips must be disposed of immediately in the household waste.

**Substances used**

Taste quality	Substance
sour	Citric acid
salty	Sodium chloride
sweet	Sucrose
bitter	Quinine hydrochloride dihydrat
Solvent	Aqua Bidest

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