ODOFIN™

Sniffin' Sticks

Threshold Test

burghardt
MESSTECHNIK

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English

Content

Description smell test with odour sticks 13
 Intended use 13
 User 13
 Age group 13
 Warnings 13
 Contraindication 14
 Patient's information for the smell test 14
 Usage of the Sniffin' Sticks 15
 Execution Threshold Test 15
 Evaluation of the results 17
 Shelf life 17
 Storage 17
 Disposal 17
 Ingredients 18
 Example 18
 Literature 19
 Evaluation threshold test (copy template) 34
 Normative data 35

**IMPORTANT
 READ CAREFULLY BEFORE USE
 KEEP FOR FUTURE USE**

Description smell test with odour sticks

ODOFIN™ „Sniffin` Sticks“ are used to investigate human olfactory performance. They have been developed in close cooperation with the association „Olfaktologie und Gustologie“ of the German Society for Ear, Nose and Throat Medicine, Head and Neck Surgery.

Intended use

The smell threshold is determined in a so called „Forced-choice-staircase procedure“. After a start concentration of the smell is found out, the dilution step is identified at which the smell can just be distinguished from non-smelling* sticks (blanks).

| Name | Font colour | Colour of stick's end | Number of sticks |
|----------------------------------|-------------|--------------------------------------|----------------------------|
| Threshold Test (n-Butanol) | Red | Red (target) Green, blue (blanks) | 16 Triplets = 48 sticks |
| Threshold Test (2-Phenylethanol) | Purple | Red (target) Green, blue (blanks) | 16 Triplets = 48 sticks |

*Filled with aqua conservans, which has its own odour

User

The test is to be performed by healthcare professionals only.

Age group

The test can be carried out starting at the age of 5 years.

Warnings

The olfactory sticks must not touch the patient's skin during the olfactory presentation. However, if the stick is inadvertently touched, it must be disposed of for hygienic reasons and because of the reduced shelf life of the odours and replaced by a new one.

All substances used in this test are not health hazardous in their used concentrations. Rarely allergic reactions show up after inhaling and skin contact. In case of allergic reaction contact the attending doctor and the manufacturer.

The user shall wear odour-free gloves during the procedure.

Use of the sticks only in an odourless or well-ventilated environment.

Only use the stick for recommended procedures and not for writing purposes.

Caution: Ask for intolerances before each test.

Store the product safely and out of the reach of children.

If the tip of the Sniffin' Sticks accidentally comes in contact with the gloves, the gloves need to be replaced to avoid contamination of other odours.

Contraindication

An allergy against one of the odorants used in the tests.

Patient's information for the smell test

The patient shall be informed about the test in the following way: Your sense of smell is to be tested. This will be done by means of odour sticks "Sniffin' Sticks". You will be blindfolded and smell the sticks which contain different odorants. The test will take approximately 15-20 minutes. All substances contained in the test are non-toxic and not harmful in the used concentrations. If you should not be able to smell and thus cannot decide on an answer, please just guess. The evaluation of the test will take into account that you have been asked to make a decision in any case.

Question: From which concentration on can you perceive an odorant? Repeatedly, three sticks will be presented to you one after another. Only one of these sticks contains the odorant, the other two contain a diluent. You have to find out the stick in which you think you have perceived the odorant.

Information for the examiner:

Patients often do not understand why they have to decide for a particular Sniffin' Stick, even though they did probably not smell anything. The patients may be given the additional explanation that this is a basic principle of modern test procedures and that accidentally correct answers are taken into account during evaluation of the test. During the evaluation it will be taken into account that they usually reach a score greater than zero, even though they did not actually smell anything, but just hit the correct answer by chance. Thus, reaching this score indicates that the patient cannot smell or cannot properly smell.

The Threshold Test may only be used in connection with the user manual, the holder and the Odofin™ "Sniffin' Sticks".

For refilling, only refill sets from the manufacturer are allowed to be used.

Usage of the Sniffin' Sticks

The examiner should always wear odourless gloves (e.g. made from cotton). When testing, the cap of only one stick should be removed. The cap has to be put back on the odour stick right after presenting the odour. Afterwards, the odour stick should immediately be placed and turned into the holder, with the tip pointing downwards.

For testing, the opened odour stick is held 2 cm centred in front of both nostrils and the patient is asked to smell, e.g. using the sentence "Please sniff now". Each stick should be presented for 3-4 seconds. The interval between the different presentations should be around 30 seconds.

If the test is to be carried out left or right sided, the opened stick is simply held in front of one nostril. The untested side is in this case closed by the patient by placing the pad of the right or left thumb on the nostril from below. The closure must not deform the nose.

To gain results as significant as possible the patients should not have anything but water for at least 15 minutes before the test. They should also refrain from consuming chewing gum, sweets or cigarettes. The actual testing should take place in a quiet and well-ventilated room. A short anamnesis can be recorded in the evaluation form. This includes: Name, age, occupation and smoking habits of the patient, as well as the information whether the test was made bilateral (L+R), left sided (L) or right sided (R). The patients should not receive any hint if their answers are correct during the test.

Execution Threshold Test

The patients are familiarised with the smell of the target stick. Therefore the stick with the highest concentration (stick no. 1 with the red stick's end) is presented. For the actual test the patient is blindfolded by putting on an "eye cover" (e.g. a sleeping mask or surgical mask). Three sticks (triplet) are each presented to the patient at approximately 5 seconds intervals. Only one stick of this triplet contains the odorant (the stick with the red cap), the other two contain just the diluent (sticks with blue and green cap = blanks). It is the patient's task to find out the stick that smells different than the two other sticks. If the target stick is recognised correctly, the triplet is presented again (see note, page 18). Each stick is presented only once, repeated

presentation of a triplet on demand of the patient is not admissible. The interval between the first stick of one triplet and the first stick of the next triplet should be 30 seconds. The highest dilution level is stick 16 (high dilution = low fragrance concentration, low dilution = high fragrance concentration). The examiner always holds three sticks in his hand while testing. In the process only the three sticks should be taken from the holder, which are currently needed. The presentation order of the stick containing the odorant has to be alternated by the examiner within each triplet. To accomplish this, the three sticks of a triplet are first presented in the order “red”, “green”, “blue”, according to the colour code of the caps. For the next triplet in the order “blue”, “red”, “green” and then in the order “green”, “blue”, “red”. This cycle is repeated throughout the whole test.

| Presentation | sequence of the sticks |
|---|------------------------|
| 1 | red – green – blue |
| 2 | blue – red – green |
| 3 | green – blue – red |
| 4 | red – green – blue |
| The sample repeats itself until the test is over! | |

1. When starting the test, triplets are presented in descending order 16, 14, 12, until the patient has recognised a triplet correctly. A certain odorant concentration **is only identified correctly if the stick containing the odorant is recognised twice in a row**, that is, if the stick containing the odorant is identified again in a second presentation of the same triplet (although a triplet is only presented a second time if the patient has identified the odorant correctly during the first presentation). **The patient always has to make a choice, even if not sure - “forced choice”**.
2. This initially correctly identified dilution step is the starting point for the following test procedure. It is marked with two crosses in the first column of the enclosed protocol sheet. After that, the next higher dilution step is presented. This is now documented in the second column.
3. If this is also identified correctly twice, again the next higher dilution step is presented and so on, until the patient makes a wrong decision. Again, the corresponding dilution step, at which the error occurred, is marked with a minus in the protocol sheet, this time in the second column.
4. After that, the next lower dilution step is presented.
5. If this is not identified by the patient, the next lower dilution step is presented and so on, until the patient identifies a dilution step correctly (i.e. twice). Again, this dilution step is marked with two crosses in the protocol sheet. Then a higher dilution step is presented and so on.
6. The test is finished when 7 turning points are passed through. The smell threshold is defined as the average value of the dilution steps which have been marked at the last 4 turning points which were passed through.

The procedure is explained again under example (see page 18).

Note: Since the blanks are not empty, but instead contain a diluent, they have a faint characteristic smell. This is intentional because the dilutions of the target sticks are made with this substance. This way a common background smell exists for target and blanks.

Evaluation of the results

The results of the Threshold Test with n-Butanol and 2-Phenylethanol show comparable results, which is shown in a comparison study ².

Since the Threshold Test by itself is not intended for a detailed assessment, the only possible statement that can be made is based on the norm values on the evaluation sheet. If the result is within the range of the norm values, there is no hint to a smell disorder of the patient.

The result of the Threshold Test is usually evaluated within the scope of the Extended Test together with the Discrimination and the Identification Test.

Shelf life

The sticks are labelled with a lot-number and a best before date. Use beyond the expiry date may lead to incorrect results.

Storage

Refer to the symbols for how to store the test (dry, protect from sunlight). The temperature should be between 10°C and 30°C.

The sticks must be stored in an upright "overhead" position (white cap downwards).

Storage in the fridge can lead to a change in the intensity of the odour.

Keep the test out of the reach of patients or children.

Disposal

Unused sticks must be disposed of in their sealed state in the household waste.

Ingredients

| | |
|----------|-------------------------------------|
| Stick | Chemical |
| No. 1-16 | 2-Phenylethanol or n-Butanol |
| No. 1-16 | Propylenglycolum or Aqua conservans |

Example

X means target detected
 - means target not detected

| Dilution | Turning points | | | | | | |
|----------|----------------|----|----|---|----|----|----|
| | 1 | 2 | 3 | 4 | 5 | 6 | 7 |
| 1 | | | | | | | |
| 2 | | | | | | | |
| 3 | | | | | | | |
| 4 | | | | | | | |
| 5 | | | | | | | |
| 6 | xx | | xx | | xx | | xx |
| 7 | ↑ | xx | - | - | | xx | - |
| 8 | x- | xx | x- | | | x- | |
| 9 | | x- | | | | | |
| 10 | - | | | | | | |
| 11 | | | | | | | |
| 12 | - | | | | | | |
| 13 | | | | | | | |
| 14 | - | | | | | | |
| 15 | | | | | | | |
| 16 | - | | | | | | |

1: In this example, the starting point is at the 6th dilution step, since the target stick was detected twice in series. Previously, steps 16, 14, 12, 10 and 8 were not detected or were not detected twice. Step 6 is also the first turning point. (Start with stick 16)

2: The next higher dilution step has now been selected (7). Two presentations were detected, i.e. the next higher dilution step (8). This was repeated until no two target sticks were detected in the 9th dilution step. This means that the second turning point has been found.

3: The next lower dilution step was then selected (8). No two target sticks were detected here, i.e. the next lower dilution level (7). No target stick was detected during the first presentation. In step 6, two target sticks have been detected. This means that the third turning point has been found.

This procedure is repeated until seven turning points have been found.

4: The score is determined from the last four turning points (i.e. turning points 1-3 are not taken into account). The average of the last four turning points is calculated for this purpose. In this example, the patient would have achieved a score of 6.75 $((7+6+8+6)/4=6.75)$.

Note: A turning point always means a change of direction.

Note: If stick 1 is also not recognised during the search for the first turning point, a smelling disorder is assumed and the score is set to 1¹.

Note: If stick 16 is correctly recognised twice in the first performance, this is the first turning point. As there are no sticks with a higher dilution level, the test is continued with stick 16.

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
Results (Average of the last four turning points)
 left _____ right _____

bilateral testing

| | | | | | | | |
|-------|--|--|--|--|--|--|--|
| Verd. | | | | | | | |
| 1 | | | | | | | |
| 2 | | | | | | | |
| 3 | | | | | | | |
| 4 | | | | | | | |
| 5 | | | | | | | |
| 6 | | | | | | | |
| 7 | | | | | | | |
| 8 | | | | | | | |
| 9 | | | | | | | |
| 10 | | | | | | | |
| 11 | | | | | | | |
| 12 | | | | | | | |
| 13 | | | | | | | |
| 14 | | | | | | | |
| 15 | | | | | | | |
| 16 | | | | | | | |

Result (Average of the last four turning points)

bilateral_____


Normative data

The following graph can be used for the evaluation. The graph was created with the data from the norm data study mentioned under "Literature". The graph is not separated by sex, as the difference of the olfactory performance between the sexes is not of a relevant magnitude for the evaluation.

The classification of normosmia and hyposmia is based on the data of the highest-performing age group (21-30) with regard to the ability to smell.

Since the test is not feasible with anosmics, 1 was set as the cut-off for anosmia¹.

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