

FACT SHEET

● LM-02 Length Measuring Unit



LM-02 Length Measuring Unit

The Length Measuring Unit LM-02 is used for the automatic measurement of the length of cylindrical objects in a range of 0-190 mm. It is possible to do single measurements or to measure a batch of objects and automatically calculate the mean and standard deviation.

Technical data

Measuring range:	0 - 190 mm
Measuring speed:	approx. 20 measurements / min
Measuring force:	adjust. from 0.4 - 2.0 N, Standard = 0.5 N
Display:	LCD TFT Display, 64.8 x 48.6 mm
Resolution:	0.02 mm
Interface:	RS 232 C, USB port
Operating voltage:	110-240 V, 50 - 60 Hz
Power consumption:	approx. 50 W
Dimensions:	600 x 300 x 140 mm (l x w x h)
Weight:	14 kg

Subject to alterations.