

## FACT SHEET

# • TBU Tobacco Blow-Out Unit



TBU Tobacco Blow-Out Unit

The Burghart TBU allows the recovery of tobacco and tubes from cigarettes.

Pressured air blows out the tobacco into a collection bag, while leaving the cigarette tube intact.

## Technical data

Cigarette Diameter:	Ø7 mm to Ø8.5 mm Format parts needed for each 0.5 mm range
Compressed Air Requirements:	Oil free & dry, adjustable pressure, min. 6 bar, approx. 100 l/min
Max. Air Pressure:	10 bar
Operation Temperature Range:	15°C to 30°C
Max. Dimensions:	224 x 120 x 29 mm
Net Weight:	363 g

Approx. 12 cigarettes per minute can be processed.