

## MUB-01 Masking Unit



MUB-01 Masking Unit

The Masking Unit MUB-01 is used to label cigarettes or filter in order to neutralize their logo or marking for blind tests in smoker panels.

An additional application is the partial to complete blocking of the ventilation zone. Labels of different dimensions may cover from 80° to 360° of the ventilation zone. To simulate the partial blocking of the ventilation through the fingers of the smoker, it is possible to use two small single labels on opposite sides. Main field of application of the MUB-01 is in the laboratory.

### Technical data

Performance:	30 cigarettes or filter per minute
Tolerance of the label position:	± 0.25 mm
Capacity of the hopper:	250 cigarettes / filter
Cigarette length:	70 - 120 mm
Cigarette diameter:	7.5 - 8.0 mm (other formats on request)
Label width:	6 - 20 mm
Label length:	6 - 30 mm
Voltage supply:	110 - 240 V, 50 - 60 Hz
Power consumption:	90 W
Compressed air:	5 bar (oil free / dry), approx. 2 l/min
Dimensions:	830 x 470 x 600 mm (L x W x H)
Weight:	40 kg

Subject to alterations.