


## RM20 CS Rebuild Service



RM20/CS (without hood) **before** rebuild



RM20/CSR (with hood) **after** rebuild

The smoking machine RM20/CS, the pneumatic unit incl. Canada Intense, the RGA-System and Puff Counter Printer are still an integral part of daily cigarette analysis in many smoking laboratories worldwide. In order to ensure the continued availability and smooth operation, we offer a rebuild service for the RM20/CS and the corresponding components. This rebuild service includes:

- every single component of the whole system will be dismantled
- all casing body parts will be re-varnished incl. a new engraving on front- and back wall
- consumables
- electrical renewal
- exchange of the mechanical memory against an electronical one
- working hours for the exchange of defective parts as well as for the overhaul or replacement of movable parts
- after completion and final assembly, the operational test according to a check list including electrical safety will be carried out.

For the Rebuild Service the RM20/CS, including components, must be send to us.